


Lipid Identification Tool

Installing and running LIT is easy as one-two-three steps. The tool provided was tested on windows and linux platforms, but for linux it is assumed you can install the required packages yourself. The installation section of this manual is meant for users that do not have anaconda installed. If you want to use linux or have anaconda installed see the `install.bat` file in the `install/` folder for the required packages. If you have any further questions, feel free to send your questions to: robbin.bouwmeester@ugent.be

Installation

1. Download the anaconda software from: <https://www.anaconda.com/download/>
2. Install anaconda by clicking on the executable you just downloaded
 - a. Make sure to add it to your path:



3. Run the `install.bat` file in the `install` folder. A black window will pop-up that will print text for a couple of minutes. When the black window disappears the required packages are installed and CALLC is ready to use.
4. To run LIT double click on `LIT.bat` in the main folder.